

# Honey Authenticity: An Integrated Approach with FlavourSpec® (GC-IMS), Visible Spectroscopy (VIS), and Molecular Biology (digital PCR)

Morozzi Pietro<sup>1</sup>, Scaramagli Sonia<sup>1</sup>, Melucci Dora<sup>2</sup>, Kassouf Nicholas<sup>2</sup>, Battaglia Fabio<sup>3</sup>, Rossini Cesare<sup>3</sup>



<sup>1</sup> COOP ITALIA Soc. Cooperativa, Via del lavoro 6/8 40033 Casalecchio di Reno, Bologna


<sup>2</sup> Department of Chemistry "Giacomo Ciamician", Alma Mater Studiorum – University of Bologna, Via Piero Gobetti 85, 40129 Bologna

<sup>3</sup> Lab Service Analytica S.R.L, Via Emilia 51/c 40011, Anzola dell'Emilia, Bologna

## 1. Objectives

- 1 Verification of Honey Authenticity: Botanical Origin
- 2 Comparison and Integration of Complementary Analytical Techniques
- 3 Identification of Anomalous Samples (Outliers)
- 4 Screening to Support Traditional Analyses with rapid and effective methods

## 2. Samples Analyzed



## 3. Integrated Multi-Instrumental Approach

a



### FlavourSpec® (GC-IMS)

#### Chemical Fingerprinting of Volatile Compounds

##### Materials and Methods

- (2.00 ± 0.02) g of raw honey in a sealed 20 mL vial, heated to 40°C for 8 minutes
- Injection of V=0.5 mL of static headspace
- Analysis time: 36 minutes per vial

- Background noise elimination (deNoise + Percentile Modules, GAS Dortmund)
- Unique pseudo-images for each honey



Results

b

### VIS Diffuse Reflectance Spectroscopy



#### Chemical Fingerprinting of Raw Honey Pigments

##### Materials and Methods

- Honey previously dissolved in a water bath, placed in a disposable plastic cuvette
- Spectral range: 360–830 nm, in integrating sphere
- Analysis time: ~ 5 seconds per cuvette

Results

- Measurement of **diffuse reflectance % vs. nm**
- **Color analysis:** parameterization and visualization of raw honey (colorimetric spaces)



c



### VIS Transmission Spectroscopy

#### Chemical Fingerprinting of Honey Pigments Diluted in Water

##### Materials and Methods

- 2 g of honey + 2 g of milliQ H<sub>2</sub>O, dissolved in a water bath, placed in a disposable plastic cuvette

- Spectral range: 360–830 nm, in transmission
- Analysis time: ~ 5 seconds per cuvette
- Measurement of **absorbance vs. nm**
- **Color analysis:** shades of orange/brown of aqueous honey extracts (colorimetric spaces)



Results

d

### Digital PCR

#### Quantification of Plant DNA

##### Materials and Methods

- **Workflow:** DNA extraction; Amplification.

- **14 molecular markers** specific to certain plant species typical of honey: 2 plastidial + 12 nuclear.
- The number of *housekeeping* gene copies is correlated to the amount of pollen in the analyzed honeys.
- In acacia honey, the percentage of species-specific DNA meets the standards.

Results

